



PUBLICATION SUMMARY

Better Accuracy, Fewer Procedures: Real-World Validation of the Mia™ Comprehensive Cervical Cancer Test

Introduction

Decades of data support HPV-based screening as the preferred strategy for cervical cancer prevention. This large-scale Norwegian study evaluated the Mel-Mont Medical Proofer 7-Type HPV mRNA E6/E7 Test, one of the core molecular components of the Mia Comprehensive Cervical Cancer Test, which also includes the cobas® HPV DNA test and the Mia by XytoTest® sample collection kit. The objective was to determine whether using mRNA testing to triage HPV DNA-positive women could improve positive predictive values (PPVs) and reduce unnecessary procedures compared with standard liquid-based cytology (Pap smears).

Key Insights

- 40% fewer unnecessary colposcopies than Pap cytology.
- Higher specificity (71% vs 48%) and 1.5× improved predictive accuracy (PPV 29% vs 19%).
- Same sensitivity (~73%) for detecting true high-grade precancer (CIN2+).
- 7 high-risk HPV types (16, 18, 31, 33, 45, 52, 58) detected—responsible for ~90% of cervical cancers worldwide.
- Demonstrates the value of combining HPV DNA and mRNA testing for better clinical clarity.

Study Design & Population:

- Participants: 58,029 women (ages 34–69) screened through Norway's national cervical program.
- HPV-positive cohort: 990 (5.6 %) tested positive by DNA assay; 962 had valid mRNA results.
- Endpoint: Histologically confirmed CIN2+ lesions.

Comparative Performance:

Decades of data support HPV-based screening as the preferred strategy for cervical cancer prevention. This large-scale Norwegian study evaluated the Mel-Mont Medical Proofer 7-Type HPV mRNA E6/E7 Test, one of the core molecular components of the Mia Comprehensive Cervical Cancer Test, which also includes the cobas® HPV DNA test and the Mia by XytoTest® sample collection kit. The objective was to determine whether using mRNA testing to triage HPV DNA-positive women could improve positive predictive values (PPVs) and reduce unnecessary procedures compared with standard liquid-based cytology (Pap smears).

Measure	Cytology (Pap)	Proofer 7-Type HPV mRNA Test	Result
Sensitivity (CIN2+)	76%	73%	Similar ability to detect true precancer.
Specificity	48%	71%	Far fewer false positives.
Positive Predictive Value	19%	29%	1.5× higher predictive accuracy.
Colposcopies per CIN2+ case	5.2%	3.1%	≈ 40% reduction in unnecessary colposcopies.

Genotype Risk:
HPV types 16, 18, 31, 33, 45, 52, and 58 were the strongest predictors of CIN2+ lesions.
Detection of E6/E7 mRNA expression provided meaningful insight into active oncogenic infections, not just the presence of HPV DNA, improving risk stratification.

Conclusion

The Mia Comprehensive Cervical Cancer Test, integrating the Proofer 7-Type HPV mRNA Test, cobas® HPV DNA Test, and Mia by XytoTest® kit, achieves the same sensitivity as cytology while providing far greater specificity. This combination enables clinicians to focus on patients truly at risk, resulting in approximately 40% fewer unnecessary colposcopies, improved diagnostic confidence, and a better patient experience.

Citation

Sørbye SW, Falang BM, Antonsen M. Performance of a 7-Type HPV mRNA Test in Triage of HPV DNA Primary Screen Positive Women Compared to Liquid-Based Cytology. *Journal of Molecular Pathology*. 2023; 4(2): 69–80. <https://doi.org/10.3390/jmp4020008>